

# IL17A Receptor Conjugated Antibody

Catalog No: #C56715

Package Size: #C56715-AF350 100ul #C56715-AF405 100ul #C56715-AF488 100ul #C56715-AF555 100ul #C56715-AF594 100ul #C56715-AF647 100ul #C56715-AF680 100ul #C56715-AF750 100ul #C56715-Biotin 100ul #C56715-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | IL17A Receptor Conjugated Antibody  |
| Host Species          | Rabbit  |
| Clonality             | Monoclonal  |
| Isotype               | Rabbit IgG  |
| Purification          | Affinity-chromatography   |
| Applications          | WB, IF  |
| Species Reactivity    | Human   |
| Specificity           | IL17A Receptor Antibody detects endogenous levels of total IL17A Receptor |
| Immunogen Description | A synthesized peptide derived from human IL17A Receptor                   |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750                    |
| Other Names           | IL-17 receptor; IL-17RA; IL17RA; interleukin-17 receptor;                 |
| Accession No.         | Uniprot:Q96F46  |
| Calculated MW         | 160kDa  |
| Storage               | Store at 4°C in dark for 6 months   |

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Receptor for IL17A and IL17F. Binds its IL17A ligand with low affinity, suggesting that additional components are involved in IL17A-induced signaling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.