

ApoH Conjugated Antibody

Catalog No: #C56716



Package Size: #C56716-AF350 100ul #C56716-AF405 100ul #C56716-AF488 100ul #C56716-AF555 100ul #C56716-AF594 100ul #C56716-AF647 100ul #C56716-AF680 100ul #C56716-AF750 100ul #C56716-Biotin 100ul #C56716-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ApoH Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Rat
Specificity	ApoH Antibody detects endogenous levels of total ApoH
Immunogen Description	A synthesized peptide derived from human ApoH
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	APC inhibitor; APOH; Apolipoprotein H; B2G1; B2GP1; B2GPI; BG;
Accession No.	Uniprot:P02749
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Binds to various kinds of negatively charged substances such as heparin, phospholipids, and dextran sulfate. May prevent activation of the intrinsic blood coagulation cascade by binding to phospholipids on the surface of damaged cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.