

ADAM15 Conjugated Antibody

Catalog No: #C56719



Package Size: #C56719-AF350 100ul #C56719-AF405 100ul #C56719-AF488 100ul #C56719-AF555 100ul #C56719-AF594 100ul #C56719-AF647 100ul #C56719-AF680 100ul #C56719-AF750 100ul #C56719-Biotin 100ul #C56719-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ADAM15 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	ADAM15 Antibody detects endogenous levels of total ADAM15
Immunogen Description	A synthesized peptide derived from human ADAM15
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Adam15; MDC15; Metargidin;
Accession No.	Uniprot:Q13444
Calculated MW	90,100kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Active metalloproteinase with gelatinolytic and collagenolytic activity. Plays a role in the wound healing process. Mediates both heterotypic intraepithelial cell/T-cell interactions and homotypic T-cell aggregation. Inhibits beta-1 integrin-mediated cell adhesion and migration of airway smooth muscle cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.