

CHRNA5 Conjugated Antibody

Catalog No: #C56723



Package Size: #C56723-AF350 100ul #C56723-AF405 100ul #C56723-AF488 100ul #C56723-AF555 100ul #C56723-AF594 100ul #C56723-AF647 100ul #C56723-AF680 100ul #C56723-AF750 100ul #C56723-Biotin 100ul #C56723-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CHRNA5 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CHRNA5 Antibody detects endogenous levels of total CHRNA5
Immunogen Description	A synthesized peptide derived from human CHRNA5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AChR; CHRNA5; LNCR2; NACHRA5;
Accession No.	Uniprot:P30532
Calculated MW	55kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.