

Tbx3 Conjugated Antibody

Catalog No: #C56724

Package Size: #C56724-AF350 100ul #C56724-AF405 100ul #C56724-AF488 100ul #C56724-AF555 100ul #C56724-AF594 100ul #C56724-AF647 100ul #C56724-AF680 100ul #C56724-AF750 100ul #C56724-Biotin 100ul #C56724-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Tbx3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Tbx3 Antibody detects endogenous levels of total Tbx3
Immunogen Description	A synthesized peptide derived from human Tbx3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL;
Accession No.	Uniprot:O15119
Calculated MW	85kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.