

WIF1 Conjugated Antibody

Catalog No: #C56726

Package Size: #C56726-AF350 100ul #C56726-AF405 100ul #C56726-AF488 100ul #C56726-AF555 100ul #C56726-AF594 100ul #C56726-AF647 100ul #C56726-AF680 100ul #C56726-AF750 100ul #C56726-Biotin 100ul #C56726-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	WIF1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	WIF1 Antibody detects endogenous levels of total WIF1
Immunogen Description	A synthesized peptide derived from human WIF1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	wif1; Wnt inhibitory factor 1;
Accession No.	Uniprot:Q9Y5W5
Calculated MW	41kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Binds to WNT proteins and inhibits their activities. May be involved in mesoderm segmentation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.