

TIP49A Conjugated Antibody

Catalog No: #C56737

Package Size: #C56737-AF350 100ul #C56737-AF405 100ul #C56737-AF488 100ul #C56737-AF555 100ul #C56737-AF594 100ul #C56737-AF647 100ul #C56737-AF680 100ul #C56737-AF750 100ul #C56737-Biotin 100ul #C56737-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TIP49A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	TIP49A Antibody detects endogenous levels of total TIP49A
Immunogen Description	A synthesized peptide derived from human TIP49A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ECP54; Pontin 52; PONTIN; RuvB like 1; RUVBL1; RVB1; TAP54 alpha; TIP49; TIP49a;
Accession No.	Uniprot:Q9Y265
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Possesses single-stranded DNA-stimulated ATPase and ATP-dependent DNA helicase (3' to 5') activity. Component of the NuA4 histone acetyltransferase complex which is involved in transcriptional activation of select genes principally by acetylation of nucleosomal histones H4 and H2A.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.