

STC1 Conjugated Antibody

Catalog No: #C56740



Package Size: #C56740-AF350 100ul #C56740-AF405 100ul #C56740-AF488 100ul #C56740-AF555 100ul #C56740-AF594 100ul

#C56740-AF647 100ul #C56740-AF680 100ul #C56740-AF750 100ul #C56740-Biotin 100ul #C56740-Conjugated 50ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	STC1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	STC1 Antibody detects endogenous levels of total STC1
Immunogen Description	A synthesized peptide derived from human STC1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Stanniocalcin; STC; STC-1; Stc1;
Accession No.	Uniprot:P52823
Calculated MW	31-37kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.