

PLCB3 Conjugated Antibody

Catalog No: #C56742

Package Size: #C56742-AF350 100ul #C56742-AF405 100ul #C56742-AF488 100ul #C56742-AF555 100ul #C56742-AF594 100ul #C56742-AF647 100ul #C56742-AF680 100ul #C56742-AF750 100ul #C56742-Biotin 100ul #C56742-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PLCB3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	PLCB3 Antibody detects endogenous levels of total PLCB3
Immunogen Description	A synthesized peptide derived from human PLCB3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Phosphoinositide phospholipase C; Phospholipase C beta 3; PLC beta 3; Plcb3;
Accession No.	Uniprot:Q01970
Calculated MW	150kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

The production of the second messenger molecules diacylglycerol (DAG) and inositol 1,4,5-trisphosphate (IP3) is mediated by activated phosphatidylinositol-specific phospholipase C enzymes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.