

LHX3 Conjugated Antibody

Catalog No: #C56743

Package Size: #C56743-AF350 100ul #C56743-AF405 100ul #C56743-AF488 100ul #C56743-AF555 100ul #C56743-AF594 100ul #C56743-AF647 100ul #C56743-AF680 100ul #C56743-AF750 100ul #C56743-Biotin 100ul #C56743-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	LHX3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	LHX3 Antibody detects endogenous levels of total LHX3
Immunogen Description	A synthesized peptide derived from human LHX3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CPHD3; LHX3; Lim3; M2 LHX3; mLim-3; mLIM3; P LIM;
Accession No.	Uniprot:Q9UBR4
Calculated MW	43kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Acts as a transcriptional activator. Binds to and activates the promoter of the alpha-glycoprotein gene, and synergistically enhances transcription from the prolactin promoter in cooperation with Pit-1.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.