

# FRA2 Conjugated Antibody

Catalog No: #C56748

Package Size: #C56748-AF350 100ul #C56748-AF405 100ul #C56748-AF488 100ul #C56748-AF555 100ul #C56748-AF594 100ul #C56748-AF647 100ul #C56748-AF680 100ul #C56748-AF750 100ul #C56748-Biotin 100ul #C56748-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FRA2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	FRA2 Antibody detects endogenous levels of total FRA2
Immunogen Description	A synthesized peptide derived from human FRA2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Fos L2; Fosl2; FRA 2;
Accession No.	Uniprot:P15408
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Controls osteoclast survival and size (By similarity). As a dimer with JUN, activates LIF transcription (By similarity). Activates CEBPB transcription in PGE2-activated osteoblasts.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.