

TMEFF2 Conjugated Antibody

Catalog No: #C56749

Package Size: #C56749-AF350 100ul #C56749-AF405 100ul #C56749-AF488 100ul #C56749-AF555 100ul #C56749-AF594 100ul #C56749-AF647 100ul #C56749-AF680 100ul #C56749-AF750 100ul #C56749-Biotin 100ul #C56749-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TMEFF2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	TMEFF2 Antibody detects endogenous levels of total TMEFF2
Immunogen Description	A synthesized peptide derived from human TMEFF2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TMEFF2; HPP1; TPEF; TR-2; Tomoregulin; TENB2; Tomoregulin-2; CT120.2;
Accession No.	Uniprot:Q9UIK5
Calculated MW	41kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

May be a survival factor for hippocampal and mesencephalic neurons. The shedded form up-regulates cancer cell proliferation, probably by promoting ERK1/2 phosphorylation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.