

# IMP2 Conjugated Antibody

Catalog No: #C56750

Package Size: #C56750-AF350 100ul #C56750-AF405 100ul #C56750-AF488 100ul #C56750-AF555 100ul #C56750-AF594 100ul #C56750-AF647 100ul #C56750-AF680 100ul #C56750-AF750 100ul #C56750-Biotin 100ul #C56750-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IMP2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	IMP2 Antibody detects endogenous levels of total IMP2
Immunogen Description	A synthesized peptide derived from human IMP2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	gf2bp2; IMP2; p62; VICKZ2;
Accession No.	Uniprot:Q9Y6M1
Calculated MW	66,62kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Binds to the 5'-UTR of the insulin-like growth factor 2 (IGF2) mRNAs. Binding is isoform-specific. May regulate translation of target mRNAs.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.