

ATP7b Conjugated Antibody

Catalog No: #C56752

Package Size: #C56752-AF350 100ul #C56752-AF405 100ul #C56752-AF488 100ul #C56752-AF555 100ul #C56752-AF594 100ul #C56752-AF647 100ul #C56752-AF680 100ul #C56752-AF750 100ul #C56752-Biotin 100ul #C56752-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ATP7b Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	ATP7b Antibody detects endogenous levels of total ATP7b
Immunogen Description	A synthesized peptide derived from human ATP7b
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ATP7B; PWD; WC1; WD; WND; WND/140 kDa;
Accession No.	Uniprot:P35670
Calculated MW	157kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Involved in the export of copper out of the cells, such as the efflux of hepatic copper into the bile.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.