

KLK2 Conjugated Antibody

Catalog No: #C56754

Package Size: #C56754-AF350 100ul #C56754-AF405 100ul #C56754-AF488 100ul #C56754-AF555 100ul #C56754-AF594 100ul #C56754-AF647 100ul #C56754-AF680 100ul #C56754-AF750 100ul #C56754-Biotin 100ul #C56754-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	KLK2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	KLK2 Antibody detects endogenous levels of total KLK2
Immunogen Description	A synthesized peptide derived from human KLK2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	hK2; hGK-1; KLK2A2;
Accession No.	Uniprot:P20151
Calculated MW	33kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Glandular kallikreins cleave Met-Lys and Arg-Ser bonds in kininogen to release Lys-bradykinin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.