

SESN1 Conjugated Antibody

Catalog No: #C56760



Package Size: #C56760-AF350 100ul #C56760-AF405 100ul #C56760-AF488 100ul #C56760-AF555 100ul #C56760-AF594 100ul #C56760-AF647 100ul #C56760-AF680 100ul #C56760-AF750 100ul #C56760-Biotin 100ul #C56760-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SESN1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	SESN1 Antibody detects endogenous levels of total SESN1
Immunogen Description	A synthesized peptide derived from human SESN1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PA26; sesn1; SEST1; sestrin 1;
Accession No.	Uniprot:Q9Y6P5
Calculated MW	70kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Involved in the reduction of peroxiredoxins. May also be regulator of cellular growth.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.