

DNAJA1 Conjugated Antibody

Catalog No: #C56764

Package Size: #C56764-AF350 100ul #C56764-AF405 100ul #C56764-AF488 100ul #C56764-AF555 100ul #C56764-AF594 100ul #C56764-AF647 100ul #C56764-AF680 100ul #C56764-AF750 100ul #C56764-Biotin 100ul #C56764-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	DNAJA1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	DNAJA1 Antibody detects endogenous levels of total DNAJA1
Immunogen Description	A synthesized peptide derived from human DNAJA1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DJ2; DjA1; DNAJ2; Dnaja1; HDJ2; HSDJ; Hsj2; HSPF4; NEDD7;
Accession No.	Uniprot:P31689
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Co-chaperone of Hsc70. Seems to play a role in protein import into mitochondria.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.