

PIM2 Conjugated Antibody

Catalog No: #C56765

Package Size: #C56765-AF350 100ul #C56765-AF405 100ul #C56765-AF488 100ul #C56765-AF555 100ul #C56765-AF594 100ul #C56765-AF647 100ul #C56765-AF680 100ul #C56765-AF750 100ul #C56765-Biotin 100ul #C56765-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PIM2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	PIM2 Antibody detects endogenous levels of total PIM2
Immunogen Description	A synthesized peptide derived from human PIM2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PIM2; Pim2h;
Accession No.	Uniprot:Q9P1W9
Calculated MW	37kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Promotes cell survival in response to a variety of proliferative signals via positive regulation of the I-kappaB kinase/NF-kappaB cascade; this process requires phosphorylation of MAP3K8/COT. Prevents apoptosis induced by growth factor withdrawal via inhibition of caspase-3 activation, and via phosphorylation of pro-apoptotic proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.