

## PR3 Conjugated Antibody

Catalog No: #C56769



Package Size: #C56769-AF350 100ul #C56769-AF405 100ul #C56769-AF488 100ul #C56769-AF555 100ul #C56769-AF594 100ul #C56769-AF647 100ul #C56769-AF680 100ul #C56769-AF750 100ul #C56769-Biotin 100ul #C56769-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PR3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	PR3 Antibody detects endogenous levels of total PR3
Immunogen Description	A synthesized peptide derived from human PR3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PRTN3; ACPA; CANCA;
Accession No.	Uniprot:P24158
Calculated MW	27kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Polymorphonuclear leukocyte serine protease that degrades elastin, fibronectin, laminin, vitronectin, and collagen types I, III, and IV (in vitro) and causes emphysema when administered by tracheal insufflation to hamsters.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.