

Cystatin A Conjugated Antibody

Catalog No: #C56770

Package Size: #C56770-AF350 100ul #C56770-AF405 100ul #C56770-AF488 100ul #C56770-AF555 100ul #C56770-AF594 100ul #C56770-AF647 100ul #C56770-AF680 100ul #C56770-AF750 100ul #C56770-Biotin 100ul #C56770-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cystatin A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Cystatin A Antibody detects endogenous levels of total Cystatin A
Immunogen Description	A synthesized peptide derived from human Cystatin A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AREI; cstA; Cystatin A; Cystatin AS; Stefin A; STF1; STFA;
Accession No.	Uniprot:P01040
Calculated MW	13kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

This is an intracellular thiol proteinase inhibitor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.