

IRF4 Conjugated Antibody

Catalog No: #C56774

Package Size: #C56774-AF350 100ul #C56774-AF405 100ul #C56774-AF488 100ul #C56774-AF555 100ul #C56774-AF594 100ul #C56774-AF647 100ul #C56774-AF680 100ul #C56774-AF750 100ul #C56774-Biotin 100ul #C56774-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IRF4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	IRF4 Antibody detects endogenous levels of total IRF4
Immunogen Description	A synthesized peptide derived from human IRF4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Interferon regulatory factor 4; Irf4; LSIRF; MUM1; NFEM5;
Accession No.	Uniprot:Q15306
Calculated MW	52kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Transcriptional activator. Binds to the interferon-stimulated response element (ISRE) of the MHC class I promoter. Binds the immunoglobulin lambda light chain enhancer, together with PU.1. Probably plays a role in ISRE-targeted signal transduction mechanisms specific to lymphoid cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.