

SENP2 Conjugated Antibody

Catalog No: #C56775



Package Size: #C56775-AF350 100ul #C56775-AF405 100ul #C56775-AF488 100ul #C56775-AF555 100ul #C56775-AF594 100ul #C56775-AF647 100ul #C56775-AF680 100ul #C56775-AF750 100ul #C56775-Biotin 100ul #C56775-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	SENP2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	SENP2 Antibody detects endogenous levels of total SENP2
Immunogen Description	A synthesized peptide derived from human SENP2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Axam2; Senp2; Sentrin specific protease 2; SMT3 specific isopeptidase 2; Smt3ip2;
Accession No.	Uniprot:Q9HC62
Calculated MW	68kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Protease that catalyzes two essential functions in the SUMO pathway: processing of full-length SUMO1, SUMO2 and SUMO3 to their mature forms and deconjugation of SUMO1, SUMO2 and SUMO3 from targeted proteins. May down-regulate CTNNB1 levels and thereby modulate the Wnt pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.