

TMOD2 Conjugated Antibody

Catalog No: #C56777

Package Size: #C56777-AF350 100ul #C56777-AF405 100ul #C56777-AF488 100ul #C56777-AF555 100ul #C56777-AF594 100ul #C56777-AF647 100ul #C56777-AF680 100ul #C56777-AF750 100ul #C56777-Biotin 100ul #C56777-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TMOD2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	TMOD2 Antibody detects endogenous levels of total TMOD2
Immunogen Description	A synthesized peptide derived from human TMOD2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Tropomodulin 2; NTMOD; TMOD2;
Accession No.	Uniprot:Q9NZR1
Calculated MW	40kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Blocks the elongation and depolymerization of the actin filaments at the pointed end. The Tmod/TM complex contributes to the formation of the short actin protofilament, which in turn defines the geometry of the membrane skeleton (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.