

CDC123 Conjugated Antibody

Catalog No: #C56778



Package Size: #C56778-AF350 100ul #C56778-AF405 100ul #C56778-AF488 100ul #C56778-AF555 100ul #C56778-AF594 100ul #C56778-AF647 100ul #C56778-AF680 100ul #C56778-AF750 100ul #C56778-Biotin 100ul #C56778-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CDC123 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CDC123 Antibody detects endogenous levels of total CDC123
Immunogen Description	A synthesized peptide derived from human CDC123
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	C10orf7; CDC123; Cell division cycle 123; D123; HT 1080; PZ32;
Accession No.	Uniprot:O75794
Calculated MW	39kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Required for S phase entry of the cell cycle.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.