

ANGPTL3 Conjugated Antibody

Catalog No: #C56782



Package Size: #C56782-AF350 100ul #C56782-AF405 100ul #C56782-AF488 100ul #C56782-AF555 100ul #C56782-AF594 100ul #C56782-AF647 100ul #C56782-AF680 100ul #C56782-AF750 100ul #C56782-Biotin 100ul #C56782-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ANGPTL3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	ANGPTL3 Antibody detects endogenous levels of total ANGPTL3
Immunogen Description	A synthesized peptide derived from human ANGPTL3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ANG5; Angiopoietin 5; Angiopoietin like 3; ANGPT5; ANGPTL3; ANL3; FHBL2;
Accession No.	Uniprot:Q9Y5C1
Calculated MW	54kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Acts in part as a hepatokine that is involved in regulation of lipid and glucose metabolism. Proposed to play a role in the trafficking of energy substrates to either storage or oxidative tissues in response to food intake (By similarity). Has a stimulatory effect on plasma triglycerides (TG), which is achieved by suppressing plasma TG clearance via inhibition of LPL activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.