

# CRMP4 Conjugated Antibody

Catalog No: #C56783

Package Size: #C56783-AF350 100ul #C56783-AF405 100ul #C56783-AF488 100ul #C56783-AF555 100ul #C56783-AF594 100ul #C56783-AF647 100ul #C56783-AF680 100ul #C56783-AF750 100ul #C56783-Biotin 100ul #C56783-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CRMP4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	CRMP4 Antibody detects endogenous levels of total CRMP4
Immunogen Description	A synthesized peptide derived from human CRMP4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CRMP4; DPYSL3; DRP3; LCRMP; TUC4; ULIP1;
Accession No.	Uniprot:Q14195
Calculated MW	62kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. Plays a role in axon guidance, neuronal growth cone collapse and cell migration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.