

PDLIM1 Conjugated Antibody

Catalog No: #C56785

Package Size: #C56785-AF350 100ul #C56785-AF405 100ul #C56785-AF488 100ul #C56785-AF555 100ul #C56785-AF594 100ul #C56785-AF647 100ul #C56785-AF680 100ul #C56785-AF750 100ul #C56785-Biotin 100ul #C56785-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PDLIM1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PDLIM1 Antibody detects endogenous levels of total PDLIM1
Immunogen Description	A synthesized peptide derived from human PDLIM1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CLIM1; CLP36; Elfin; Enigma; hCLIM1; Pdlim1;
Accession No.	Uniprot:O00151
Calculated MW	36kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Cytoskeletal protein that may act as an adapter that brings other proteins (like kinases) to the cytoskeleton.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.