

# HIRA Conjugated Antibody

Catalog No: #C56787

Package Size: #C56787-AF350 100ul #C56787-AF405 100ul #C56787-AF488 100ul #C56787-AF555 100ul #C56787-AF594 100ul #C56787-AF647 100ul #C56787-AF680 100ul #C56787-AF750 100ul #C56787-Biotin 100ul #C56787-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HIRA Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	HIRA Antibody detects endogenous levels of total HIRA
Immunogen Description	A synthesized peptide derived from human HIRA
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DGCR1; HIR; HIRA; TUP 1; TUPLE1;
Accession No.	Uniprot:P54198
Calculated MW	112kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Cooperates with ASF1A to promote replication-independent chromatin assembly. Required for the periodic repression of histone gene transcription during the cell cycle.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.