

CHX10 Conjugated Antibody

Catalog No: #C56791

Package Size: #C56791-AF350 100ul #C56791-AF405 100ul #C56791-AF488 100ul #C56791-AF555 100ul #C56791-AF594 100ul #C56791-AF647 100ul #C56791-AF680 100ul #C56791-AF750 100ul #C56791-Biotin 100ul #C56791-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CHX10 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CHX10 Antibody detects endogenous levels of total CHX10
Immunogen Description	A synthesized peptide derived from human CHX10
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CHX 10; HOX10; MCOP2; MCOPCB3; RET1; Vsx2;
Accession No.	Uniprot:P58304
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Plays a significant role in the specification and morphogenesis of the sensory retina. May also participate in the development of the cells of the inner nuclear layer, particularly bipolar cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.