

FXR1 Conjugated Antibody

Catalog No: #C56795



Package Size: #C56795-AF350 100ul #C56795-AF405 100ul #C56795-AF488 100ul #C56795-AF555 100ul #C56795-AF594 100ul #C56795-AF647 100ul #C56795-AF680 100ul #C56795-AF750 100ul #C56795-Biotin 100ul #C56795-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	FXR1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	FXR1 Antibody detects endogenous levels of total FXR1
Immunogen Description	A synthesized peptide derived from human FXR1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	FXR1; FXR1P; hFXR1p;
Accession No.	Uniprot:P51114
Calculated MW	70-85kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

RNA-binding protein required for embryonic and postnatal development of muscle tissue. May regulate intracellular transport and local translation of certain mRNAs.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.