

SSTR5 Conjugated Antibody

Catalog No: #C56796



Package Size: #C56796-AF350 100ul #C56796-AF405 100ul #C56796-AF488 100ul #C56796-AF555 100ul #C56796-AF594 100ul

#C56796-AF647 100ul #C56796-AF680 100ul #C56796-AF750 100ul #C56796-Biotin 100ul #C56796-Conjugated 50ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	SSTR5 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	SSTR5 Antibody detects endogenous levels of total SSTR5
Immunogen Description	A synthesized peptide derived from human SSTR5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	SSTR5; SS5R; Sstr5;
Accession No.	Uniprot:P35346
Calculated MW	42kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Receptor for somatostatin 28 and to a lesser extent for somatostatin-14. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.