

FKBP38 Conjugated Antibody

Catalog No: #C56798

Package Size: #C56798-AF350 100ul #C56798-AF405 100ul #C56798-AF488 100ul #C56798-AF555 100ul #C56798-AF594 100ul #C56798-AF647 100ul #C56798-AF680 100ul #C56798-AF750 100ul #C56798-Biotin 100ul #C56798-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FKBP38 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	FKBP38 Antibody detects endogenous levels of total FKBP38
Immunogen Description	A synthesized peptide derived from human FKBP38
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	FKBP 38; Fkbp8; FKBP38; hFKBP38; Rotamase; Sam11;
Accession No.	Uniprot:Q14318
Calculated MW	52kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Constitutively inactive PPIase, which becomes active when bound to calmodulin and calcium. Seems to act as a chaperone for BCL2, targets it to the mitochondria and modulates its phosphorylation state.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.