

CELF1 Conjugated Antibody

Catalog No: #C56800

Package Size: #C56800-AF350 100ul #C56800-AF405 100ul #C56800-AF488 100ul #C56800-AF555 100ul #C56800-AF594 100ul #C56800-AF647 100ul #C56800-AF680 100ul #C56800-AF750 100ul #C56800-Biotin 100ul #C56800-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CELF1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	CELF1 Antibody detects endogenous levels of total CELF1
Immunogen Description	A synthesized peptide derived from human CELF1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BRUNOL2; celf1; CUGBP; CUGBP1; EDEN BP; hNab50; NAB50; NAPOR;
Accession No.	Uniprot:Q92879
Calculated MW	52kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

RNA-binding protein implicated in the regulation of several post-transcriptional events. Involved in pre-mRNA alternative splicing, mRNA translation and stability. Mediates exon inclusion and/or exclusion in pre-mRNA that are subject to tissue-specific and developmentally regulated alternative splicing.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.