

# TCF4 Conjugated Antibody

Catalog No: #C56801

Package Size: #C56801-AF350 100ul #C56801-AF405 100ul #C56801-AF488 100ul #C56801-AF555 100ul #C56801-AF594 100ul #C56801-AF647 100ul #C56801-AF680 100ul #C56801-AF750 100ul #C56801-Biotin 100ul #C56801-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TCF4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	TCF4 Antibody detects endogenous levels of total TCF4
Immunogen Description	A synthesized peptide derived from human TCF4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	E2-2; ITF2; PTHS; SEF2; ITF-2; SEF-2; TCF-4; SEF2-1; SEF2-1A; SEF2-1B; bHLHb19;
Accession No.	Uniprot:P15884
Calculated MW	71,90kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Transcription factor that binds to the immunoglobulin enhancer Mu-E5/KE5-motif. Binds to the E-box present in the somatostatin receptor 2 initiator element (SSTR2-INR) to activate transcription (By similarity). Preferentially binds to either 5'-ACANNTGT-3' or 5'-CCANNTGG-3'.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.