

# RhoGAP Conjugated Antibody

Catalog No: #C56804

Package Size: #C56804-AF350 100ul #C56804-AF405 100ul #C56804-AF488 100ul #C56804-AF555 100ul #C56804-AF594 100ul #C56804-AF647 100ul #C56804-AF680 100ul #C56804-AF750 100ul #C56804-Biotin 100ul #C56804-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | RhoGAP Conjugated Antibody  |
| Host Species          | Rabbit  |
| Clonality             | Monoclonal  |
| Isotype               | Rabbit IgG  |
| Purification          | Affinity-chromatography   |
| Applications          | WB, IF, FC  |
| Species Reactivity    | Human Mouse Rat   |
| Specificity           | RhoGAP Antibody detects endogenous levels of total RhoGAP                       |
| Immunogen Description | A synthesized peptide derived from human RhoGAP                                 |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750                          |
| Other Names           | ARHGAP5; p190 B; p190 RhoGAP; Rho GTPase activating protein 5; RhoGAP; RhoGAP5; |
| Accession No.         | Uniprot:Q13017  |
| Calculated MW         | 190kDa  |
| Storage               | Store at 4°C in dark for 6 months   |

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

GTPase-activating protein for Rho family members (PubMed:8537347).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.