

PI 3 Kinase R4 Conjugated Antibody

Catalog No: #C56805

Package Size: #C56805-AF350 100ul #C56805-AF405 100ul #C56805-AF488 100ul #C56805-AF555 100ul #C56805-AF594 100ul #C56805-AF647 100ul #C56805-AF680 100ul #C56805-AF750 100ul #C56805-Biotin 100ul #C56805-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PI 3 Kinase R4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PI 3 Kinase R4 Antibody detects endogenous levels of total PI 3 Kinase R4
Immunogen Description	A synthesized peptide derived from human PI 3 Kinase R4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	P150; Phosphoinositide 3 kinase regulatory subunit 4; PI3 kinase p150 subunit; PIK3R4; VPS15;
Accession No.	Uniprot:Q99570
Calculated MW	150kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Regulatory subunit of the PI3K complex. May regulate membrane trafficking late in the endocytic pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.