

## eIF3e Conjugated Antibody

Catalog No: #C56806



Package Size: #C56806-AF350 100ul #C56806-AF405 100ul #C56806-AF488 100ul #C56806-AF555 100ul #C56806-AF594 100ul #C56806-AF647 100ul #C56806-AF680 100ul #C56806-AF750 100ul #C56806-Biotin 100ul #C56806-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | eIF3e Conjugated Antibody                               |
| Host Species          | Rabbit  |
| Clonality             | Monoclonal  |
| Isotype               | Rabbit IgG  |
| Purification          | Affinity-chromatography                                 |
| Applications          | WB, IF  |
| Species Reactivity    | Human Mouse Rat   |
| Specificity           | eIF3e Antibody detects endogenous levels of total eIF3e |
| Immunogen Description | A synthesized peptide derived from human eIF3e          |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750  |
| Other Names           | eIF3e; EIF3S6; eIFe;                                    |
| Accession No.         | Uniprot:P60228  |
| Calculated MW         | 45kDa   |
| Storage               | Store at 4°C in dark for 6 months                       |

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Component of the eukaryotic translation initiation factor 3 (eIF-3) complex, which is required for several steps in the initiation of protein synthesis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.