

# KDM1 / LSD1 Conjugated Antibody

Catalog No: #C56807

Package Size: #C56807-AF350 100ul #C56807-AF405 100ul #C56807-AF488 100ul #C56807-AF555 100ul #C56807-AF594 100ul #C56807-AF647 100ul #C56807-AF680 100ul #C56807-AF750 100ul #C56807-Biotin 100ul #C56807-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	KDM1 / LSD1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	KDM1 / LSD1 Antibody detects endogenous levels of total KDM1 / LSD1
Immunogen Description	A synthesized peptide derived from human KDM1 / LSD1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AOF2; CPRF; EC1; KDM1; Kdm1a; LSD1;
Accession No.	Uniprot:O60341
Calculated MW	110kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Histone demethylase that demethylates both 'Lys-4' (H3K4me) and 'Lys-9' (H3K9me) of histone H3, thereby acting as a coactivator or a corepressor, depending on the context.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.