

Wnt5b Conjugated Antibody

Catalog No: #C56808

Package Size: #C56808-AF350 100ul #C56808-AF405 100ul #C56808-AF488 100ul #C56808-AF555 100ul #C56808-AF594 100ul #C56808-AF647 100ul #C56808-AF680 100ul #C56808-AF750 100ul #C56808-Biotin 100ul #C56808-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Wnt5b Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Wnt5b Antibody detects endogenous levels of total Wnt5b
Immunogen Description	A synthesized peptide derived from human Wnt5b
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Wnt5b;
Accession No.	Uniprot:Q9H1J7
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.