

# SLC22A3 Conjugated Antibody

Catalog No: #C56809

Package Size: #C56809-AF350 100ul #C56809-AF405 100ul #C56809-AF488 100ul #C56809-AF555 100ul #C56809-AF594 100ul #C56809-AF647 100ul #C56809-AF680 100ul #C56809-AF750 100ul #C56809-Biotin 100ul #C56809-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	SLC22A3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	SLC22A3 Antibody detects endogenous levels of total SLC22A3
Immunogen Description	A synthesized peptide derived from human SLC22A3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	EMT; EMTH; OCT3; Orct3; Slc22a3;
Accession No.	Uniprot:O75751
Calculated MW	75kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Mediates potential-dependent transport of a variety of organic cations. May play a significant role in the disposition of cationic neurotoxins and neurotransmitters in the brain.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.