

CXCL5 Conjugated Antibody

Catalog No: #C56810

Package Size: #C56810-AF350 100ul #C56810-AF405 100ul #C56810-AF488 100ul #C56810-AF555 100ul #C56810-AF594 100ul #C56810-AF647 100ul #C56810-AF680 100ul #C56810-AF750 100ul #C56810-Biotin 100ul #C56810-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CXCL5 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	CXCL5 Antibody detects endogenous levels of total CXCL5
Immunogen Description	A synthesized peptide derived from human CXCL5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AMCFII; Cxcl5; ENA78; SCYB5;
Accession No.	Uniprot:P42830
Calculated MW	10kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Involved in neutrophil activation. In vitro, ENA-78(8-78) and ENA-78(9-78) show a threefold higher chemotactic activity for neutrophil granulocytes.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.