

ARF5 Conjugated Antibody

Catalog No: #C56812



Package Size: #C56812-AF350 100ul #C56812-AF405 100ul #C56812-AF488 100ul #C56812-AF555 100ul #C56812-AF594 100ul #C56812-AF647 100ul #C56812-AF680 100ul #C56812-AF750 100ul #C56812-Biotin 100ul #C56812-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ARF5 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	ARF5 Antibody detects endogenous levels of total ARF5
Immunogen Description	A synthesized peptide derived from human ARF5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ADP ribosylation factor 5; Arf5;
Accession No.	Uniprot:P84085
Calculated MW	18kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

GTP-binding protein that functions as an allosteric activator of the cholera toxin catalytic subunit, an ADP-ribosyltransferase. Involved in protein trafficking; may modulate vesicle budding and uncoating within the Golgi apparatus.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.