

ABI2 Conjugated Antibody

Catalog No: #C56814

Package Size: #C56814-AF350 100ul #C56814-AF405 100ul #C56814-AF488 100ul #C56814-AF555 100ul #C56814-AF594 100ul #C56814-AF647 100ul #C56814-AF680 100ul #C56814-AF750 100ul #C56814-Biotin 100ul #C56814-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ABI2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	ABI2 Antibody detects endogenous levels of total ABI2
Immunogen Description	A synthesized peptide derived from human ABI2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Abelson interactor 2; ABI2; ABI2B; Abl binding protein 3; AblBP3; ArgBP1; ARGBPIA; ArgBPIB; SSH3BP2;
Accession No.	Uniprot:Q9NYB9
Calculated MW	75kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May act in regulation of cell growth and transformation by interacting with nonreceptor tyrosine kinases ABL1 and/or ABL2. Part of the WAVE complex that regulates lamellipodia formation. The WAVE complex regulates actin filament reorganization via its interaction with the Arp2/3 complex. Regulates ABL1/c-Abl-mediated phosphorylation of MENA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.