

HLA F Conjugated Antibody

Catalog No: #C56816

Package Size: #C56816-AF350 100ul #C56816-AF405 100ul #C56816-AF488 100ul #C56816-AF555 100ul #C56816-AF594 100ul #C56816-AF647 100ul #C56816-AF680 100ul #C56816-AF750 100ul #C56816-Biotin 100ul #C56816-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	HLA F Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	HLA F Antibody detects endogenous levels of total HLA F
Immunogen Description	A synthesized peptide derived from human HLA F
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HLAF;
Accession No.	Uniprot:P30511
Calculated MW	42kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Involved in the presentation of foreign antigens to the immune system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.