

Neurocan Conjugated Antibody

Catalog No: #C56817



Package Size: #C56817-AF350 100ul #C56817-AF405 100ul #C56817-AF488 100ul #C56817-AF555 100ul #C56817-AF594 100ul

#C56817-AF647 100ul #C56817-AF680 100ul #C56817-AF750 100ul #C56817-Biotin 100ul #C56817-Conjugated 50ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Neurocan Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Neurocan Antibody detects endogenous levels of total Neurocan
Immunogen Description	A synthesized peptide derived from human Neurocan
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Cspg3; Ncan; NEUR; Neurocan core protein;
Accession No.	Uniprot:O14594
Calculated MW	130,270kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May modulate neuronal adhesion and neurite growth during development by binding to neural cell adhesion molecules (NG-CAM and N-CAM).
Chondroitin sulfate proteoglycan; binds to hyaluronic acid.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.