

# Eph receptor B3 Conjugated Antibody

Catalog No: #C56827

Package Size: #C56827-AF350 100ul #C56827-AF405 100ul #C56827-AF488 100ul #C56827-AF555 100ul #C56827-AF594 100ul #C56827-AF647 100ul #C56827-AF680 100ul #C56827-AF750 100ul #C56827-Biotin 100ul #C56827-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Eph receptor B3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Eph receptor B3 Antibody detects endogenous levels of total Eph receptor B3
Immunogen Description	A synthesized peptide derived from human Eph receptor B3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Cek10; EK2; Embryonic kinase 2; ephb3; ETK2; hEK2; Mdk5; Sek4; TYRO6;
Accession No.	Uniprot:P54753
Calculated MW	110kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Receptor for members of the ephrin-B family. Binds to ephrin-B1 and -B2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.