

# HMGCS2 Conjugated Antibody

Catalog No: #C56833

Package Size: #C56833-AF350 100ul #C56833-AF405 100ul #C56833-AF488 100ul #C56833-AF555 100ul #C56833-AF594 100ul #C56833-AF647 100ul #C56833-AF680 100ul #C56833-AF750 100ul #C56833-Biotin 100ul #C56833-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HMGCS2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	HMGCS2 Antibody detects endogenous levels of total HMGCS2
Immunogen Description	A synthesized peptide derived from human HMGCS2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HMG CoA synthase; HMGCS2;
Accession No.	Uniprot:P54868
Calculated MW	50kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

This enzyme condenses acetyl-CoA with acetoacetyl-CoA to form HMG-CoA, which is the substrate for HMG-CoA reductase.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.