

# VPS28 Conjugated Antibody

Catalog No: #C56841

Package Size: #C56841-AF350 100ul #C56841-AF405 100ul #C56841-AF488 100ul #C56841-AF555 100ul #C56841-AF594 100ul #C56841-AF647 100ul #C56841-AF680 100ul #C56841-AF750 100ul #C56841-Biotin 100ul #C56841-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	VPS28 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	VPS28 Antibody detects endogenous levels of total VPS28
Immunogen Description	A synthesized peptide derived from human VPS28
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	H-Vps28; Vps28;
Accession No.	Uniprot:Q9UK41
Calculated MW	25kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Component of the ESCRT-I complex, a regulator of vesicular trafficking process.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.