

VEGFD Conjugated Antibody

Catalog No: #C56842

Package Size: #C56842-AF350 100ul #C56842-AF405 100ul #C56842-AF488 100ul #C56842-AF555 100ul #C56842-AF594 100ul #C56842-AF647 100ul #C56842-AF680 100ul #C56842-AF750 100ul #C56842-Biotin 100ul #C56842-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	VEGFD Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	VEGFD Antibody detects endogenous levels of total VEGFD
Immunogen Description	A synthesized peptide derived from human VEGFD
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	FIGF; VEGFD;
Accession No.	Uniprot:O43915
Calculated MW	40kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Growth factor active in angiogenesis, lymphangiogenesis and endothelial cell growth, stimulating their proliferation and migration and also has effects on the permeability of blood vessels.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.