

KLC1 Conjugated Antibody

Catalog No: #C56849

Package Size: #C56849-AF350 100ul #C56849-AF405 100ul #C56849-AF488 100ul #C56849-AF555 100ul #C56849-AF594 100ul #C56849-AF647 100ul #C56849-AF680 100ul #C56849-AF750 100ul #C56849-Biotin 100ul #C56849-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	KLC1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	KLC1 Antibody detects endogenous levels of total KLC1
Immunogen Description	A synthesized peptide derived from human KLC1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Kinesin2; KLC; Klc1; KNS2; KNS2A;
Accession No.	Uniprot:Q07866
Calculated MW	65kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Kinesin is a microtubule-associated force-producing protein that may play a role in organelle transport. The light chain may function in coupling of cargo to the heavy chain or in the modulation of its ATPase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.